

Antibodies Product Datasheet

Product Name:	ATF-6 β Antibody
Catalog Number:	C11585
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATF6B. AA range:401-450
Concentration:	1 mg/ml
Gene Name:	ATF6B
Protein Name:	Cyclic AMP-dependent transcription factor ATF-6 beta
Other Name:	ATF6B; CREBL1; G13; Cyclic AMP-dependent transcription factor ATF-6 beta; cAMP-dependent transcription factor ATF-6 beta; Activating transcription factor 6 beta; ATF6-beta; Protein G13; cAMP response element-binding protein-related protein;
Human Swiss Prot No.:	Q99941
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ATF-6 β Polyclonal Antibody detects endogenous levels of ATF-6 β protein.
Storage and Stability:	-20°C/1 year