

Antibodies Product Datasheet

Product Name:	MMTAG2 Antibody
Catalog Number:	C11613
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MMTAG2. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	MMTAG2
Protein Name:	Multiple myeloma tumor-associated protein 2
Other Name:	MMTAG2; C1orf35; Multiple myeloma tumor-associated protein 2; hMMTAG2
Human Swiss Prot No.:	Q9BU76
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MMTAG2 Polyclonal Antibody detects endogenous levels of MMTAG2 protein.
Storage and Stability:	-20°C/1 year