

Antibodies Product Datasheet

Product Name:	Kanadaptin Antibody
Catalog Number:	C11628
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NADAP. AA range:421-470
Concentration:	1 mg/ml
Gene Name:	SLC4A1AP
Protein Name:	Kanadaptin
Other Name:	SLC4A1AP; HLC3; Kanadaptin; Human lung cancer oncogene 3 protein; HLC-3; Kidney anion exchanger adapter protein; Solute carrier family 4 anion exchanger member 1 adapter protein
Human Swiss Prot No.:	Q9BWU0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Kanadaptin Polyclonal Antibody detects endogenous levels of Kanadaptin protein.
Storage and Stability:	-20°C/1 year