

Antibodies Product Datasheet

Product Name:	SPZ1 Antibody
Catalog Number:	C11631
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SPZ1. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	SPZ1
Protein Name:	Spermatogenic leucine zipper protein 1
Other Name:	SPZ1; TSP1; Spermatogenic leucine zipper protein 1; Testis-specific protein 1; Testis-specific protein NYD-TSP1
Human Swiss Prot No.:	Q9BXG8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SPZ1 Polyclonal Antibody detects endogenous levels of SPZ1 protein.
Storage and Stability:	-20°C/1 year