

Antibodies Product Datasheet

Product Name:	DEC2 Antibody
Catalog Number:	C11661
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BHLHB3. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	BHLHE41
Protein Name:	Class E basic helix-loop-helix protein 41
Other Name:	BHLHE41; BHLHB3; DEC2; SHARP1; Class E basic helix-loop-helix protein 41; bHLHe41; Class B basic helix-loop-helix protein 3; bHLHb3; Differentially expressed in chondrocytes protein 2; hDEC2; Enhancer-of-split and hairy-related protein 1; S
Human Swiss Prot No.:	Q9C0J9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DEC2 Polyclonal Antibody detects endogenous levels of 40879 protein.
Storage and Stability:	-20°C/1 year