

Antibodies Product Datasheet

Product Name:	PEG3 Antibody
Catalog Number:	C11667
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PEG3. AA range:1031-1080
Concentration:	1 mg/ml
Gene Name:	PEG3
Protein Name:	Paternally-expressed gene 3 protein
Other Name:	PEG3; KIAA0287; ZSCAN24; Paternally-expressed gene 3 protein; Zinc finger and SCAN domain-containing protein 24
Human Swiss Prot No.:	Q9GZU2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PEG3 Polyclonal Antibody detects endogenous levels of PEG3 protein.
Storage and Stability:	-20°C/1 year