

Antibodies Product Datasheet

Product Name:	ECRG4 Antibody
Catalog Number:	C11679
Reactivity:	Human, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ECRG4. AA range:41-90
Concentration:	1 mg/ml
Gene Name:	C2orf40
Protein Name:	Augurin
Other Name:	C2orf40; ECRG4; Augurin; Esophageal cancer-related gene 4 protein
Human Swiss Prot No.:	Q9H1Z8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ECRG4 Polyclonal Antibody detects endogenous levels of ECRG4 protein.
Storage and Stability:	-20°C/1 year