

Antibodies Product Datasheet

Product Name:	Eps8L2 Antibody
Catalog Number:	C11699
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ES8L2. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	EPS8L2
Protein Name:	Epidermal growth factor receptor kinase substrate 8-like protein 2
Other Name:	EPS8L2; EPS8R2; PP13181; Epidermal growth factor receptor kinase substrate 8-like protein 2; EPS8-like protein 2; Epidermal growth factor receptor pathway substrate 8-related protein 2; EPS8-related protein 2
Human Swiss Prot No.:	Q9H6S3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Eps8L2 Polyclonal Antibody detects endogenous levels of Eps8L2 protein.
Storage and Stability:	-20°C/1 year