

Antibodies Product Datasheet

Product Name:	InsP5 2-kinase Antibody
Catalog Number:	C11708
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IPPK. AA range:11-60
Concentration:	1 mg/ml
Gene Name:	IPPK
Protein Name:	Inositol-pentakisphosphate 2-kinase
Other Name:	IPPK; C9orf12; Inositol-pentakisphosphate 2-kinase; IPK1 homolog; Inositol-1; 3,4,5,6-pentakisphosphate 2-kinase; Ins(1,3,4,5,6)P5 2-kinase; InsP5 2-kinase
Human Swiss Prot No.:	Q9H8X2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	InsP5 2-kinase Polyclonal Antibody detects endogenous levels of InsP5 2-kinase protein.
Storage and Stability:	-20°C/1 year