

Antibodies Product Datasheet

Product Name:	Arnt 2 Antibody
Catalog Number:	C11723
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ARNT2. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	ARNT2
Protein Name:	Aryl hydrocarbon receptor nuclear translocator 2
Other Name:	ARNT2; BHLHE1; KIAA0307; Aryl hydrocarbon receptor nuclear translocator 2; ARNT protein 2; Class E basic helix-loop-helix protein 1; bHLHe1
Human Swiss Prot No.:	Q9HBZ2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Arnt 2 Polyclonal Antibody detects endogenous levels of Arnt 2 protein.
Storage and Stability:	-20°C/1 year