

Antibodies Product Datasheet

Product Name:	BCA3 Antibody
Catalog Number:	C11741
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BCA3. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	AKIP1
Protein Name:	A-kinase-interacting protein 1
Other Name:	AKIP1; BCA3; C11orf17; A-kinase-interacting protein 1; Breast cancer-associated gene 3 protein; PKA-interacting protein; Proline-rich protein BCA3
Human Swiss Prot No.:	Q9NQ31
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BCA3 Polyclonal Antibody detects endogenous levels of BCA3 protein.
Storage and Stability:	-20°C/1 year