

Antibodies Product Datasheet

Product Name:	CHRAC15 Antibody
Catalog Number:	C11750
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CHRC1. AA range:81-130
Concentration:	1 mg/ml
Gene Name:	CHRAC1
Protein Name:	Chromatin accessibility complex protein 1
Other Name:	CHRAC1; CHRAC15; Chromatin accessibility complex protein 1; CHRAC-1; Chromatin accessibility complex 15 kDa protein; CHRAC-15; HuCHRAC15; DNA polymerase epsilon subunit p15
Human Swiss Prot No.:	Q9NRG0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CHRAC15 Polyclonal Antibody detects endogenous levels of CHRAC15 protein.
Storage and Stability:	-20°C/1 year