

Antibodies Product Datasheet

Product Name:	Kpm Antibody
Catalog Number:	C11751
Reactivity:	Human, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LATS2. AA range:541-590
Concentration:	1 mg/ml
Gene Name:	LATS2
Protein Name:	Serine/threonine-protein kinase LATS2
Other Name:	LATS2; KPM; Serine/threonine-protein kinase LATS2; Kinase phosphorylated during mitosis protein; Large tumor suppressor homolog 2; Serine/threonine-protein kinase kpm; Warts-like kinase
Human Swiss Prot No.:	Q9NRM7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Kpm Polyclonal Antibody detects endogenous levels of Kpm protein.
Storage and Stability:	-20°C/1 year