

Antibodies Product Datasheet

Product Name:	CRSP77 Antibody
Catalog Number:	C11762
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MED17. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	MED17
Protein Name:	Mediator of RNA polymerase II transcription subunit 17
Other Name:	MED17; ARC77; CRSP6; DRIP77; DRIP80; TRAP80; Mediator of RNA polymerase II transcription subunit 17; Activator-recruited cofactor 77 kDa component; ARC77; Cofactor required for Sp1 transcriptional activation subunit 6; CRSP complex subunit
Human Swiss Prot No.:	Q9NVC6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CRSP77 Polyclonal Antibody detects endogenous levels of CRSP77 protein.
Storage and Stability:	-20°C/1 year