

## Antibodies Product Datasheet

|                               |                                                                                                                       |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | Ethanolamine kinase 2 Antibody                                                                                        |
| <b>Catalog Number:</b>        | C11763                                                                                                                |
| <b>Reactivity:</b>            | Human, Mouse                                                                                                          |
| <b>Applications:</b>          | WB, IHC-p, IF(paraffin section), ELISA                                                                                |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.                                    |
| <b>Modification:</b>          | Unmodified/Total                                                                                                      |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                     |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Immunogen:</b>             | The antiserum was produced against synthesized peptide derived from human EK12. AA range:51-100                       |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                               |
| <b>Gene Name:</b>             | ETNK2                                                                                                                 |
| <b>Protein Name:</b>          | Ethanolamine kinase 2                                                                                                 |
| <b>Other Name:</b>            | ETNK2; EK12; HMFT1716; Ethanolamine kinase 2; EK1 2; Ethanolamine kinase-like protein                                 |
| <b>Human Swiss Prot No.:</b>  | Q9NVF9                                                                                                                |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| <b>Specificity:</b>           | Ethanolamine kinase 2 Polyclonal Antibody detects endogenous levels of Ethanolamine kinase 2 protein.                 |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                          |