

Antibodies Product Datasheet

Product Name:	ASF1B Antibody
Catalog Number:	C11765
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ASF1B. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	ASF1B
Protein Name:	Histone chaperone ASF1B
Other Name:	ASF1B; Histone chaperone ASF1B; Anti-silencing function protein 1 homolog B; hAsf1; hAsf1b; CCG1-interacting factor A-II; CIA-II; hCIA-II
Human Swiss Prot No.:	Q9NVP2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ASF1B Polyclonal Antibody detects endogenous levels of ASF1B protein.
Storage and Stability:	-20°C/1 year