

Antibodies Product Datasheet

Product Name:	DUS2L Antibody
Catalog Number:	C11772
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DUS2L. AA range:421-470
Concentration:	1 mg/ml
Gene Name:	DUS2L
Protein Name:	tRNA-dihydrouridine(20) synthase [NAD(P)+]-like
Other Name:	DUS2L; DUS2; tRNA-dihydrouridine(20) synthase [NAD(P)+]-like; Up-regulated in lung cancer protein 8; URLC8; tRNA-dihydrouridine synthase 2-like; hDUS2
Human Swiss Prot No.:	Q9NX74
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DUS2L Polyclonal Antibody detects endogenous levels of DUS2L protein.
Storage and Stability:	-20°C/1 year