

Antibodies Product Datasheet

Product Name:	MLTK Antibody
Catalog Number:	C11781
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MLTK. AA range:701-750
Concentration:	1 mg/ml
Gene Name:	MLTK
Protein Name:	Mitogen-activated protein kinase kinase kinase MLT
Other Name:	MLTK; ZAK; HCCS4; Mitogen-activated protein kinase kinase kinase MLT; Human cervical cancer suppressor gene 4 protein; HCCS-4; Leucine zipper- and sterile alpha motif-containing kinase; MLK-like mitogen-activated protein triple kinase; Mixe
Human Swiss Prot No.:	Q9NYL2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MLTK Polyclonal Antibody detects endogenous levels of MLTK protein.
Storage and Stability:	-20°C/1 year