

Antibodies Product Datasheet

Product Name:	FASTKD2 Antibody
Catalog Number:	C11785
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FAKD2. AA range:171-220
Concentration:	1 mg/ml
Gene Name:	FASTKD2
Protein Name:	FAST kinase domain-containing protein 2
Other Name:	FASTKD2; KIAA0971; FAST kinase domain-containing protein 2
Human Swiss Prot No.:	Q9NYY8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FASTKD2 Polyclonal Antibody detects endogenous levels of FASTKD2 protein.
Storage and Stability:	-20°C/1 year