

Antibodies Product Datasheet

Product Name:	BRAF35 Antibody
Catalog Number:	C11798
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HMG20B. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	HMG20B
Protein Name:	SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily E member 1-related
Other Name:	HMG20B; BRAF35; HMGX2; HMGXB2; SMARCE1R; SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily E member 1-related; SMARCE1-related protein; BRCA2-associated factor 35; HMG box-containing protein 20B; HMG domain-Q9P0W2
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BRAF35 Polyclonal Antibody detects endogenous levels of BRAF35 protein.
Storage and Stability:	-20°C/1 year