

Antibodies Product Datasheet

Product Name:	Dnmt3b Antibody
Catalog Number:	C11809
Reactivity:	Human, Chicken(testedbyyourcustomer)
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DNMT3B. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	DNMT3B
Protein Name:	DNA (cytosine-5)-methyltransferase 3B
Other Name:	DNMT3B; DNA (cytosine-5)-methyltransferase 3B; Dnmt3b; DNA methyltransferase HsaIIIB; DNA MTase HsaIIIB; M.HsaIIIB
Human Swiss Prot No.:	Q9UBC3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Dnmt3b Polyclonal Antibody detects endogenous levels of Dnmt3b protein.
Storage and Stability:	-20°C/1 year