

Antibodies Product Datasheet

Product Name:	Dbf4 Antibody
Catalog Number:	C11821
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DBF4. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	DBF4
Protein Name:	Protein DBF4 homolog A
Other Name:	DBF4; ASK; DBF4A; ZDBF1; Protein DBF4 homolog A; Activator of S phase kinase; Chiffon homolog A; DBF4-type zinc finger-containing protein 1
Human Swiss Prot No.:	Q9UBU7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Dbf4 Polyclonal Antibody detects endogenous levels of Dbf4 protein.
Storage and Stability:	-20°C/1 year