

Antibodies Product Datasheet

Product Name:	MRG15 Antibody
Catalog Number:	C11822
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MORF4L1. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	MORF4L1
Protein Name:	Mortality factor 4-like protein 1
Other Name:	MORF4L1; MRG15; FWP006; HSPC008; HSPC061; PP368; Mortality factor 4-like protein 1; MORF-related gene 15 protein; Protein MSL3-1; Transcription factor-like protein MRG15
Human Swiss Prot No.:	Q9UBU8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MRG15 Polyclonal Antibody detects endogenous levels of MRG15 protein.
Storage and Stability:	-20°C/1 year