

Antibodies Product Datasheet

Product Name:	EphA6 Antibody
Catalog Number:	C11829
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EPHA6. AA range:891-940
Concentration:	1 mg/ml
Gene Name:	EPHA6
Protein Name:	Ephrin type-A receptor 6
Other Name:	EPHA6; EHK2; HEK12; Ephrin type-A receptor 6; EPH homology kinase 2; EHK-2; EPH-like kinase 12; EK12
Human Swiss Prot No.:	Q9UF33
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EphA6 Polyclonal Antibody detects endogenous levels of EphA6 protein.
Storage and Stability:	-20°C/1 year