

Antibodies Product Datasheet

Product Name:	CGGBP1 Antibody
Catalog Number:	C11831
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CGGBP1. AA range:41-90
Concentration:	1 mg/ml
Gene Name:	CGGBP1
Protein Name:	CGG triplet repeat-binding protein 1
Other Name:	CGGBP1; CGGBP; CGG triplet repeat-binding protein 1; CGG-binding protein 1; 20 kDa CGG-binding protein; p20-CGGBP DNA-binding protein
Human Swiss Prot No.:	Q9UFW8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CGGBP1 Polyclonal Antibody detects endogenous levels of CGGBP1 protein.
Storage and Stability:	-20°C/1 year