

Antibodies Product Datasheet

Product Name:	ARMCX3 Antibody
Catalog Number:	C11839
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ARMX3. AA range:291-340
Concentration:	1 mg/ml
Gene Name:	ARMCX3
Protein Name:	Armadillo repeat-containing X-linked protein 3
Other Name:	ARMCX3; ALEX3; BM-017; Armadillo repeat-containing X-linked protein 3; ARM protein lost in epithelial cancers on chromosome X 3; Protein ALEX3
Human Swiss Prot No.:	Q9UH62
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ARMCX3 Polyclonal Antibody detects endogenous levels of ARMCX3 protein.
Storage and Stability:	-20°C/1 year