

Antibodies Product Datasheet

Product Name:	TRPS1 Antibody
Catalog Number:	C11843
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TRPS1. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	TRPS1
Protein Name:	Zinc finger transcription factor Trps1
Other Name:	TRPS1; Zinc finger transcription factor Trps1; Tricho-rhino-phalangeal syndrome type I protein; Zinc finger protein GC79
Human Swiss Prot No.:	Q9UHF7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TRPS1 Polyclonal Antibody detects endogenous levels of TRPS1 protein.
Storage and Stability:	-20°C/1 year