

Antibodies Product Datasheet

Product Name:	IP6K2 Antibody
Catalog Number:	C11845
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IP6K2. AA range:161-210
Concentration:	1 mg/ml
Gene Name:	IP6K2
Protein Name:	Inositol hexakisphosphate kinase 2
Other Name:	IP6K2; IHPK2; TCCCA00113; Inositol hexakisphosphate kinase 2; InsP6 kinase 2; P(i)-uptake stimulator; PiUS
Human Swiss Prot No.:	Q9UHH9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IP6K2 Polyclonal Antibody detects endogenous levels of IP6K2 protein.
Storage and Stability:	-20°C/1 year