

Antibodies Product Datasheet

Product Name:	Ikaros 3 Antibody
Catalog Number:	C11873
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IKZF3. AA range:361-410
Concentration:	1 mg/ml
Gene Name:	IKZF3
Protein Name:	Zinc finger protein Aiolos
Other Name:	IKZF3; ZNFN1A3; Zinc finger protein Aiolos; Ikaros family zinc finger protein 3
Human Swiss Prot No.:	Q9UKT9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Ikaros 3 Polyclonal Antibody detects endogenous levels of Ikaros 3 protein.
Storage and Stability:	-20°C/1 year