

Antibodies Product Datasheet

Product Name:	Endomucin Antibody
Catalog Number:	C11879
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Mucin-14. AA range:201-250
Concentration:	1 mg/ml
Gene Name:	EMCN
Protein Name:	Endomucin
Other Name:	EMCN; EMCN2; MUC14; Endomucin; Endomucin-2; Gastric cancer antigen Ga34; Mucin-14; MUC-14
Human Swiss Prot No.:	Q9ULC0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Endomucin Polyclonal Antibody detects endogenous levels of Endomucin protein.
Storage and Stability:	-20°C/1 year