

Antibodies Product Datasheet

Product Name:	CRSP130 Antibody
Catalog Number:	C11881
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MED23. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	MED23
Protein Name:	Mediator of RNA polymerase II transcription subunit 23
Other Name:	MED23; ARC130; CRSP3; DRIP130; KIAA1216; SUR2; Mediator of RNA polymerase II transcription subunit 23; Activator-recruited cofactor 130 kDa component; ARC130; Cofactor required for Sp1 transcriptional activation subunit 3; CRSP complex subu
Human Swiss Prot No.:	Q9ULK4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CRSP130 Polyclonal Antibody detects endogenous levels of CRSP130 protein.
Storage and Stability:	-20°C/1 year