

Antibodies Product Datasheet

Product Name:	JIP-3 Antibody
Catalog Number:	C11903
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human JIP3. AA range:621-670
Concentration:	1 mg/ml
Gene Name:	MAPK8IP3
Protein Name:	C-Jun-amino-terminal kinase-interacting protein 3
Other Name:	MAPK8IP3; JIP3; KIAA1066; C-Jun-amino-terminal kinase-interacting protein 3; JIP-3; JNK-interacting protein 3; JNK MAP kinase scaffold protein 3; Mitogen-activated protein kinase 8-interacting protein 3
Human Swiss Prot No.:	Q9UPT6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	JIP-3 Polyclonal Antibody detects endogenous levels of JIP-3 protein.
Storage and Stability:	-20°C/1 year