

Antibodies Product Datasheet

Product Name:	DLEC1 Antibody
Catalog Number:	C11919
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DLEC1. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	DLEC1
Protein Name:	Deleted in lung and esophageal cancer protein 1
Other Name:	DLEC1; DLC1; Deleted in lung and esophageal cancer protein 1; Deleted in lung cancer protein 1; DLC-1
Human Swiss Prot No.:	Q9Y238
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DLEC1 Polyclonal Antibody detects endogenous levels of DLEC1 protein.
Storage and Stability:	-20°C/1 year