

Antibodies Product Datasheet

Product Name:	ChoKB Antibody
Catalog Number:	C11923
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CHKB. AA range:301-350
Concentration:	1 mg/ml
Gene Name:	CHKB
Protein Name:	Choline/ethanolamine kinase
Other Name:	CHKB; CHETK; CHKL; Choline/ethanolamine kinase; Choline kinase beta; CK; CKB; Choline kinase-like protein; Ethanolamine kinase; EK; Ethanolamine kinase beta; EKB; choline/ethanolamine kinase beta; CKEKB
Human Swiss Prot No.:	Q9Y259
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ChoKB Polyclonal Antibody detects endogenous levels of ChoKB protein.
Storage and Stability:	-20°C/1 year