

Antibodies Product Datasheet

Product Name:	Asf1a Antibody
Catalog Number:	C11928
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ASF1A. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	ASF1A
Protein Name:	Histone chaperone ASF1A
Other Name:	ASF1A; CGI-98; HSPC146; Histone chaperone ASF1A; Anti-silencing function protein 1 homolog A; hAsf1; hAsf1a; CCG1-interacting factor A; CIA; hCIA
Human Swiss Prot No.:	Q9Y294
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Asf1a Polyclonal Antibody detects endogenous levels of Asf1a protein.
Storage and Stability:	-20°C/1 year