

Antibodies Product Datasheet

Product Name:	ATF-5 Antibody
Catalog Number:	C11931
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATF5. AA range:221-270
Concentration:	1 mg/ml
Gene Name:	ATF5
Protein Name:	Cyclic AMP-dependent transcription factor ATF-5
Other Name:	ATF5; ATFX; Cyclic AMP-dependent transcription factor ATF-5; cAMP-dependent transcription factor ATF-5; Activating transcription factor 5; Transcription factor ATFX
Human Swiss Prot No.:	Q9Y2D1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ATF-5 Polyclonal Antibody detects endogenous levels of ATF-5 protein.
Storage and Stability:	-20°C/1 year