

Antibodies Product Datasheet

Product Name:	AK6 Antibody
Catalog Number:	C11948
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KAD6. AA range:11-60
Concentration:	1 mg/ml
Gene Name:	
Protein Name:	Adenylate kinase isoenzyme 6
Other Name:	TAF9; AK6; CINAP; AD-004; CGI-137; Adenylate kinase isoenzyme 6; ATP-AMP transphosphorylase 6; Coilin-interacting nuclear ATPase protein; hCINAP
Human Swiss Prot No.:	Q9Y3D8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AK6 Polyclonal Antibody detects endogenous levels of AK6 protein.
Storage and Stability:	-20°C/1 year