

Antibodies Product Datasheet

Product Name:	Dyrk1B Antibody
Catalog Number:	C11955
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DYR1B. AA range:331-380
Concentration:	1 mg/ml
Gene Name:	DYRK1B
Protein Name:	Dual specificity tyrosine-phosphorylation-regulated kinase 1B
Other Name:	DYRK1B; MIRK; Dual specificity tyrosine-phosphorylation-regulated kinase 1B; Minibrain-related kinase; Mirk protein kinase
Human Swiss Prot No.:	Q9Y463
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Dyrk1B Polyclonal Antibody detects endogenous levels of Dyrk1B protein.
Storage and Stability:	-20°C/1 year