

Antibodies Product Datasheet

Product Name:	Matriptase Antibody
Catalog Number:	C11979
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ST14. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	ST14
Protein Name:	Suppressor of tumorigenicity 14 protein
Other Name:	ST14; PRSS14; SNC19; TADG15; Suppressor of tumorigenicity 14 protein; Matriptase; Membrane-type serine protease 1; MT-SP1; Prostamin; Serine protease 14; Serine protease TADG-15; Tumor-associated differentially-expressed gene 15 protein
Human Swiss Prot No.:	Q9Y5Y6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Matriptase Polyclonal Antibody detects endogenous levels of Matriptase protein.
Storage and Stability:	-20°C/1 year