

Antibodies Product Datasheet

Product Name:	IRAK-M Antibody
Catalog Number:	C11980
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IRAK3. AA range:491-540
Concentration:	1 mg/ml
Gene Name:	IRAK3
Protein Name:	Interleukin-1 receptor-associated kinase 3
Other Name:	IRAK3; Interleukin-1 receptor-associated kinase 3; IRAK-3; IL-1 receptor-associated kinase M; IRAK-M
Human Swiss Prot No.:	Q9Y616
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IRAK-M Polyclonal Antibody detects endogenous levels of IRAK-M protein.
Storage and Stability:	-20°C/1 year