

Antibodies Product Datasheet

Product Name:	SMRTE Antibody
Catalog Number:	C11981
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NCOR2. AA range:511-560
Concentration:	1 mg/ml
Gene Name:	NCOR2
Protein Name:	Nuclear receptor corepressor 2
Other Name:	NCOR2; CTG26; Nuclear receptor corepressor 2; N-CoR2; CTG repeat protein 26; SMAP270; Silencing mediator of retinoic acid and thyroid hormone receptor; SMRT; T3 receptor-associating factor; TRAC; Thyroid-; retinoic-acid-receptor-associated
Human Swiss Prot No.:	Q9Y618
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SMRTE Polyclonal Antibody detects endogenous levels of SMRTE protein.
Storage and Stability:	-20°C/1 year