

Antibodies Product Datasheet

Product Name:	EID-1 Antibody
Catalog Number:	C11984
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EID1. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	EID1
Protein Name:	EP300-interacting inhibitor of differentiation 1
Other Name:	EID1; C15orf3; CRI1; RBP21; PNAS-22; PTD014; EP300-interacting inhibitor of differentiation 1; 21 kDa pRb-associated protein; CREBBP/EP300 inhibitory protein 1; E1A-like inhibitor of differentiation 1; EID-1
Human Swiss Prot No.:	Q9Y6B2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EID-1 Polyclonal Antibody detects endogenous levels of EID-1 protein.
Storage and Stability:	-20°C/1 year