

Antibodies Product Datasheet

Product Name:	MEK Kinase-4 Antibody
Catalog Number:	C11995
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MAP3K4. AA range:1281-1330
Concentration:	1 mg/ml
Gene Name:	MAP3K4
Protein Name:	Mitogen-activated protein kinase kinase kinase 4
Other Name:	MAP3K4; KIAA0213; MAPKKK4; MEKK4; MTK1; Mitogen-activated protein kinase kinase kinase 4; MAP three kinase 1; MAPK/ERK kinase kinase 4; MEK kinase 4; MEKK 4
Human Swiss Prot No.:	Q9Y6R4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MEK Kinase-4 Polyclonal Antibody detects endogenous levels of MEK Kinase-4 protein.
Storage and Stability:	-20°C/1 year