

Antibodies Product Datasheet

Product Name:	SR-1F Antibody
Catalog Number:	C12012
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human 5-HT-1F. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	HTR1F
Protein Name:	5-hydroxytryptamine receptor 1F
Other Name:	HTR1F; HTR1EL; 5-hydroxytryptamine receptor 1F; 5-HT-1F; 5-HT1F; Serotonin receptor 1F
Human Swiss Prot No.:	P30939
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SR-1F Polyclonal Antibody detects endogenous levels of SR-1F protein.
Storage and Stability:	-20°C/1 year