

Antibodies Product Datasheet

Product Name:	A Cyclase IV Antibody
Catalog Number:	C12034
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ADCY4. AA range:195-244
Concentration:	1 mg/ml
Gene Name:	ADCY4
Protein Name:	Adenylate cyclase type 4
Other Name:	ADCY4; Adenylate cyclase type 4; ATP pyrophosphate-lyase 4; Adenylate cyclase type IV; Adenylyl cyclase 4
Human Swiss Prot No.:	Q8NFM4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	A Cyclase IV Polyclonal Antibody detects endogenous levels of A Cyclase IV protein.
Storage and Stability:	-20°C/1 year