

Antibodies Product Datasheet

Product Name:	A Cyclase V/VI Antibody
Catalog Number:	C12035
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ADCY5/6. AA range:931-980
Concentration:	1 mg/ml
Gene Name:	ADCY5/ADCY6
Protein Name:	Adenylate cyclase type 5/6
Other Name:	ADCY5; Adenylate cyclase type 5; ATP pyrophosphate-lyase 5; Adenylate cyclase type V; Adenylyl cyclase 5; ADCY6; KIAA0422; Adenylate cyclase type 6; ATP pyrophosphate-lyase 6; Adenylate cyclase type VI; Adenylyl cyclase 6; Ca(2+)-inhibitabl
Human Swiss Prot No.:	O95622/O43306
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	A Cyclase V/VI Polyclonal Antibody detects endogenous levels of A Cyclase V/VI protein.
Storage and Stability:	-20°C/1 year