

Antibodies Product Datasheet

Product Name:	A Cyclase VII Antibody
Catalog Number:	C12037
Reactivity:	Human, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ADCY7. AA range:191-240
Concentration:	1 mg/ml
Gene Name:	ADCY7
Protein Name:	Adenylate cyclase type 7
Other Name:	ADCY7; KIAA0037; Adenylate cyclase type 7; ATP pyrophosphate-lyase 7; Adenylate cyclase type VII; Adenylyl cyclase 7
Human Swiss Prot No.:	P51828
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	A Cyclase VII Polyclonal Antibody detects endogenous levels of A Cyclase VII protein.
Storage and Stability:	-20°C/1 year