

Antibodies Product Datasheet

Product Name:	AHA-1 Antibody
Catalog Number:	C12046
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AHSA1. AA range:161-210
Concentration:	1 mg/ml
Gene Name:	AHSA1
Protein Name:	Activator of 90 kDa heat shock protein ATPase homolog 1
Other Name:	AHSA1; C14orf3; HSPC322; Activator of 90 kDa heat shock protein ATPase homolog 1; AHA1; p38
Human Swiss Prot No.:	O95433
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AHA-1 Polyclonal Antibody detects endogenous levels of AHA-1 protein.
Storage and Stability:	-20°C/1 year