

Antibodies Product Datasheet

Product Name:	AMID Antibody
Catalog Number:	C12047
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AIFM2. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	AIFM2
Protein Name:	Apoptosis-inducing factor 2
Other Name:	AIFM2; AMID; PRG3; Apoptosis-inducing factor 2; Apoptosis-inducing factor homologous mitochondrion-associated inducer of death; Apoptosis-inducing factor-like mitochondrion-associated inducer of death; p53-responsive gene 3 protein
Human Swiss Prot No.:	Q9BRQ8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AMID Polyclonal Antibody detects endogenous levels of AMID protein.
Storage and Stability:	-20°C/1 year