

Antibodies Product Datasheet

Product Name:	API5 Antibody
Catalog Number:	C12052
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human API-5. AA range:421-470
Concentration:	1 mg/ml
Gene Name:	API5
Protein Name:	Apoptosis inhibitor 5
Other Name:	API5; MIG8; Apoptosis inhibitor 5; API-5; Antiapoptosis clone 11 protein; AAC-11; Cell migration-inducing gene 8 protein; Fibroblast growth factor 2-interacting factor; FIF; Protein XAGL
Human Swiss Prot No.:	Q9BZZ5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	API5 Polyclonal Antibody detects endogenous levels of API5 protein.
Storage and Stability:	-20°C/1 year